Registered ISO 9001:2008

Registered ISO 14001:2004

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BIOTECHNOLOGY SCIENCE

www.centrion.co.kr

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Control System	Built-in SCADA System	<ul style="list-style-type: none"> * Built-in SCADA System - Voltage Specifications : 90~260V 50/60Hz Free Voltage - Built-in type SMPS Module UL certification. - CE product certification to enable component selection and PCB ARTWORK. - PWM frequency control AC Servo Motor, Ac induction Motor, BLDC Motor and Is, slow down as fast and smooth implementation. - RS232 x 2, RS422 or RS485 communication port, USB port, apply - Analog Input : 12 points, analog output : 12 points, analog record:12 points - Temp. PH, DO Cascade, ORP, OD, O2, Co2, Agitation, MFC, Pressure Balancer Control
	Feed Control Mode	Fed-Batch Culture by DO, pH Interlock pump Control 4 x Built-in Feeding Pump(Boxer or Watson-Malow) External pump 2ea use
Communication port	Ethernet	1PC 1-6 connected to one controller (after completion of standard controller progress) Data logging, trend graph PC Control: process control (PID, the upper and lower values , p rograms, cascade, Feed)
	Record Output	Each sensor can be output by selecting data. D-SUB 25Pin Female Type, USB Excel file stored separately.
	USB	Measured data, setup data stored in the USB
Peristaltic pump	Pump	4built-in pumps, two external pump (optional)
	Motor type	AC Motor or DC Motor, minimum speed is 1rpm
	Speed range	0~70rpm
	Resolution	1rpm
	Control mode	Programmable PID Feeding control. Pump can be assigned for Acid, Base, Antifoam, Feed
Bottom Heating Plate	Range	Up to 90
	Resolution	0.1
	Power source	100-120V~50/60Hz or 210-230V~50/60Hz With electrical safety cutoff switch



Ordering information

Cat. NO.	Description
FSBC-LF-L01	LAB Fermentor Controller

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Vessel	Total Volume(Liter)	1.5L	3L	5L	7L	10L	15L
	Working Volume(Liter)	1	2.2	3	5	7	12
	Material	Borosilicate glass Autoclavable, SUS316L Stainless steel for Topplate and all fittings					
	Etc.	Single & Double 1.5L~15L Total Vol. pH, DO, Foam,Level,Pressure Probe, Addition Ports, Exhaust Port, Baffles 316L,Condenser					
Aeration	Flow rate	0~10LPM			0~20LPM		
	Option	Air Flow Meter / Mass Flow Controller / Mass Flow Manual					
	Sparger	Standard : Ring Sparger / Micro Sparger					
	Inlet Filter	0.2 μ m Disposable Hydrophobic Filter					
Agitation	Drive	Direct Top Drive Servo Motor 200W~400W,BLDC Motor / Single Mechanical seal					
	Range	10~1500rpm					
	Impellers	Rushton Standard With Fermentation / Pitched Blade Standard With Cell Culture Marine Blade or Spin Filter Note : Customized impellers are available					
Temperature	Thermostat system	0~150°C \pm 0.1°C / pt100, Ω Probe Heating&Cooling PID Control / Built-in Heat Exchanger / Automatic Cooling Water Valve					
pH	Range / Resolution	2~14.0 pH of set point, PID control / Gel Type					
	Probe	Electrode Autoclavable, InPro 3030 – MettlerToledo or Hamilton (temperature range – 0 ~ 140° C, the maximum pressure – 6bar)					
	Control mode	Programmable PID Control system					
DO	Range / Resolution	0~200% or 0~20ppm / 0.1%					
	Probe	Electrode Autoclavable Polarographic or Galvanic Oxygen Sensor – Mettler Toledo or Hamilton					
	Control mode	Programmable PID Control system DO Cascade to Agitation, Mas Flow, Feeding Pump Control Oxygen enrichment module (optional) Gas Mixing Station module (optional)					
ORP	Range / Resolution	Measuring range – 1000 ~ –1000 mv. / 1mV					
	Probe	Electrode Autoclavable Redox potential, Oxidation of measurements–Mettler Toledo or Hamilton (Temperature range – 0 ~ 140° C, the maximum pressure – 2.5bar)					
Anti Foam	Range / Sensor	Conductivity 0 ~ 300k Ω (Measuring the amount of foam)					
OD(Optional)	Range / Sensor	Measuring range 0 ... 100 EBC 0 ... 100 EBC0 ... 400 FTU (Turbidity measurement)					
MFM or MFC(Optional)	Range / Sensor	0~10L/min (Air flow measurement) / Mass Flow Meter / Mass Flow Control					
Level Sensor(Optional)	Control	Electrode type Hi/Low Vessel Level control					
Balancer(Optional)	Control	9,999.99g / RS232C Measured by the amount of weight on the scale output					
Analog Input & output	Control	Various sensors are available by utilizing the analog input / output port					
Analog Recording	Control	It can be the sensor output data value via the analog recording					

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


- Durable stainless steel supporting rods and bottom plate are designed.
- Single glass type : Pyrex, Stainless316L.
- Fast cooling/heating is achieved through the inner cooling coil and heating base unit.

Bowl Vessel	Cat. NO.	Total Vol. (Working Vol.)	Inner Dia. (mm) D	Inner Height (mm) H	Electrode Length (mm)	H:D
	FSBC-LF03-04	7L(5L)	Ø160	330	325	1.8:1
	FSBC-LF03-05	10L(7L)	Ø190	355	325	1.8:1
	FSBC-LF03-06	15L(12L)	Ø190	500	425	2.6:1

- Stainless Double Jacket in upper body, Single Glass in under body.
- Circulating water in double jacket for effective temperature control by large contact surface.
- Connected to extra Water Bath for temperature.

Double Vessel	Cat. NO.	Total Vol. (Working Vol.)	Inner Dia. (mm) D	Inner Height (mm) H	Electrode Length (mm)	H:D
	FSBC-LF02-02	3L(2.5L)	Ø130	220	255	1.7:1
	FSBC-LF02-03	5L(3L)	Ø145	300	225	1.3:1
	FSBC-LF02-04	7L(5L)	Ø160	330	325	1.8:1
	FSBC-LF02-05	10L(7L)	Ø190	355	325	1.8:1
	FSBC-LF02-06	15L(12L)	Ø190	500	425	2.6:1

- Double glass jacketed type vessel is specially designed for temperature sensitive.
- Circulating water in double jacket for effective temperature control by large contact surface.
- Connected to extra Water Bath for temperature control.
- Provides sophisticated temperature control.

Single Round Vessel	Cat. NO.	Total Vol. (Working Vol.)	Inner Dia. (mm) D	Inner Height (mm) H	Electrode Length (mm)	H:D
	FSBC-LF04-01	1.5L(1L)	Ø115	135	120	1.2:1
	FSBC-LF04-02	3L(2.5L)	Ø133	220	225	1.7:1
	FSBC-LF04-03	5L(3L)	Ø170	220	225	1.3:1
	FSBC-LF04-04	7L(5L)	Ø170	310	325	1.8:1
	FSBC-LF04-05	10L(7L)	Ø190	355	325	1.8:1

- Vessel with an inner coil to achieve a fast cooling.
- Single round vessel: round type under body.
- Temperature control: Glass surrounding heating blanket.
- Usually applied in animal cell fermentation.

PBR Vessel (Photo type)	Cat. NO.	Total Vol. (Working Vol.)	Inner Dia. (mm) D	Inner Height (mm) H	Electrode Length (mm)	H:D
	FSBC-LF05-02	3L(2.5L)	Ø133	220	225	1.7:1
	FSBC-LF05-03	5L(3L)	Ø170	220	225	1.3:1
	FSBC-LF05-04	7L(5L)	Ø170	310	325	1.8:1
	FSBC-LF05-05	10L(7L)	Ø190	355	325	1.8:1
	FSBC-LF05-06	15L(10L)	Ø190	355	325	1.8:1

- Single vessel type, LED Bar inside to make growth by exposing light into algae.
- Illumination : 430nm, 680nm, White, Red, Blue LED.
- LED inside / fluorescent jacket.

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Part Name	Cat. NO.	Description
Rushton 6-Blade Impeller	FSBC-LF-IM01-01	Rushton 6-Blade Impeller Ø50, 1,5L iter Vessel, 3ea/pk
	FSBC-LF-IM01-02	Rushton 6-Blade Impeller Ø60, 3L iter Vessel, 3ea/pk
	FSBC-LF-IM01-03	Rushton 6-Blade Impeller Ø72, 5~7L iter Vessel, 3ea/pk
	FSBC-LF-IM01-04	Rushton 6-Blade Impeller Ø80, 10~15L iter Vessel, 3ea/pk

- Rushton turbine impeller is located from middle to bottom and basically used to agitate liquid in vessel.
- It makes the flow to the same direction with turning radius and has features of paddle and propeller types to be available for high-speed rotation by minimizing resistance.
- General model to be used due to the highest ventilation efficiency.
- Available for the agitator assemblies, both lipseal.

Part Name	Cat. NO.	Description
Marine Impeller	FSBC-LF-IM04-01	Marine Impeller Ø60, 3Liter Vessel, 1ea/pk
	FSBC-LF-IM04-02	Marine Impeller Ø80, 5~7Liter Vessel, 1ea/pk
	FSBC-LF-IM04-03	Marine Impeller Ø100, 10~15Liter Vessel, 1ea/pk

- Applied to damageable cell.
- It stirs liquid from down to up rotating in low-speed since it has the large surface of wings.

Part Name	Cat. NO.	Description
Pitched Blade Impeller	FSBC-LF-IM02-01	Pitched Blade Impeller Ø60, 3Liter Vessel, 2ea/pk
	FSBC-LF-IM02-02	Pitched Blade Impeller Ø72, 5~7Liter Vessel, 2ea/pk
	FSBC-LF-IM02-03	Pitched Blade Impeller Ø80, 10~15Liter Vessel, 2ea/pk

- More efficient than flat blade impeller.
- Flow is discharged both axially and radially depending on the angle.

Part Name	Cat. NO.	Description
Hollowed Paddle Impeller	FSBC-LF-IM05-01	Hollowed Paddle Impeller Ø60, 3Liter Vessel, 1ea/pk
	FSBC-LF-IM05-02	Hollowed Paddle Marine Impeller Ø80, 5~7Liter Vessel, 1ea/pk
	FSBC-LF-IM05-03	Hollowed Paddle Marine Impeller Ø100, 10~15Liter Vessel, 1ea/pk

- It stirs liquid from up to down rotating in low-speed since it has the large surface of wings.
- Extra additional it.

Part Name	Cat. NO.	Description
Pitched Paddle Impeller	FSBC-LF-IM03-01	Pitched Paddle Impeller Ø65, 3Liter Vessel, 2ea/pk
	FSBC-LF-IM03-02	Pitched Paddle Impeller Ø85, 5~7Liter Vessel, 2ea/pk
	FSBC-LF-IM03-03	Pitched Paddle Impeller Ø10, 10~15Liter Vessel, 2ea/pk

- More efficient than flat blade impeller.
- Flow is discharged both axially and radially depending on the angle.

Part Name	Cat. NO.	Description
Foam Breaker	FSBC-LF-IM06-01	Foam Breaker 70, 3Liter Vessel, 1ea/pk
	FSBC-LF-IM06-02	Foam Breaker 90, 5~7Liter Vessel, 1ea/pk
	FSBC-LF-IM06-03	Foam Breaker 110, 10~15Liter Vessel, 1ea/pk

- Our customized foam breaker helps you to maintain the proper foam level.
- Installed in upper side of inside shaft, used in bursting fermentation foam.
- It is not used in the case of animal cell fermentation since there is barely any foam.

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Specification


Light intensity control range	10 ~ 100%
Light module	Fluorescent lamp or LED lamp
Light Color	White, RED, Blue, Green
Max. Lux	Max. 10000lux
Sensor Control & Recording	pH, DO, Temp, Foam, Pressure etc.
Sparger & Gas Mode	Micro sparger, 2 gas and 4 gas mixer(optional)
Filter & Condenser	0.2µm air filter cartridge, SUS 316
Temp Control System	Cartridge heater or chiller
Feeding Control System	Peristaltic pump control Acid, Base, Foam, Feed
Light Control Mode	On/off or set of the lighting controller
Power Voltage	110V~220V, 50/60Hz, 10A

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Cat. NO.	Description
FSBC-LF-PB01	Photo Bioreactor Autoclavable Fermentor 3L ~ 15L

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
Part Name	Cat. NO.	Description
pH Sensor	FSBC-BS-PS01	Measuring range pH 0.00 – 14.00
		<ul style="list-style-type: none"> • Mettler Toledo or Hamilton • Autoclavable /130°C for 30min • Electrode Length120~425(mm)

Part Name	Cat. NO.	Description
DO Sensor	FSBC-BS-DS01	Measuring range 0.0-20.0ppm or 0.0-200%
		<ul style="list-style-type: none"> • Mettler Toledo or Hamilton • Autoclavable /130°C for 30min • Electrode Length120~425(mm)

Part Name	Cat. NO.	Description
Temp Sensor	FSBC-BS-TS01	Measuring range 0 – 150°C
		<ul style="list-style-type: none"> • High accuracy Pt100 Sensor • Electrode Length150~450(mm)

Part Name	Cat. NO.	Description
Foam Sensor	FSBC-BS-FS01	Measuring range Conductivity 0 ~ 300k Ω (Measuring the amount of foam)
		<ul style="list-style-type: none"> • Adjustable Height • Auto/Manual/Stop Control • Stainless Steel Trips and Teflon Body • Autoclavable /130°C for 30min

Part Name	Cat. NO.	Description
Gas Mixer	FSBC-BS-GM01	Control Gas : Air, N ₂ , CO ₂ , O ₂ Gas Mixing Control parameter : DO and pH
		<ul style="list-style-type: none"> • Adjustable Height • Auto/Manual Control • 4Solenoid valves, 4Pressure gauges and rotameters

Part Name	Cat. NO.	Description
MFC/MFM	FSBC-BS-MS01	Measuring range 0 – 10L/min
		<ul style="list-style-type: none"> • Mass Flow Control (Auto control) • Digital display for easy confirmation of control • Within 0.5sec. At setting $\pm 5\%$F.Selectrode • 24VDC$\pm 10\%$ / 250mA or Less • BIOCANVAS LF Connect Auto Control

Part Name	Cat. NO.	Description
Gas Sensor	FSBC-BS-GS01	Measuring range 0.1 – 25 Vol.% O ₂ , , 1 – 50 Vol.% O ₂ , 0 – 25 Vol.% CO ₂
		<ul style="list-style-type: none"> • parallel measurement of O₂ and CO₂ • Compact stainless steel housing • PAT conform in-situ-measurement • Maker by BlueInOne • BIOCANVAS LF Connect Auto Control

Part Name	Cat. NO.	Description
Balancer	FSBC-BS-BS01	Measuring range 210g~1100g
		<ul style="list-style-type: none"> • Digital display for easy confirmation of control • The best of magnetic force Restoration • Single point parallelogram load Sensing • RS-232C Interface • BIOCANVAS LF Connect Auto Control

Part Name	Cat. NO.	Description
External Pump	FSBC-BS-EP01	Speed 0.1rpm-100rpm, CW/CCW
		<ul style="list-style-type: none"> • Digital display for easy confirmation of control • Applicable pump heads • Membrane keypad, easy to operate • BIOCANVAS LF Connect Auto Control

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Vessel volume	Total Capacity(Liter)	30L	50L	70L	100L	150L	300L	500L	1000L
	Working Volume(Liter)	21L	35L	50L	70L	100L	210L	350L	700L
Agitation	Drive & Range	30L ~500L : AC motor Top drive(mechanical seal driving) 50~1,000 ±1 RPM(In the condition of no load)							
	Control	Manual or automatic control of agitation speed							
	Impellers	Rushton turbine impeller, Foam breaker(SUS 316L)							
	Sensor & Control	Magnetic hall sensor, PID control							
	Agitation Type	Top Drive(Mechanical Seal) or Bottom Drive(Magnetic Coupling)							
Head Lift	The head lift system is powered by air cylinder and guided by two sets of linear motion guide								
Materials & Finish	Vessel/Piping	SUS 316L, SUS 304(electrolytic polishing)							
Sterilization	Heat exchanger or inject steam into jacket.(105°C~130°C)								
Aeration	Sparger	Ring sparger(SUS 316L)							
	Filter	Inlet(0,2µm), Outlet filter(0,2µm)							
pH	Range / Electrode	2~14,0 pH of set point, PID control / Gel Type							
DO	Range / Electrode	0~200% or 0~20ppm, PID control / Polarographic Type							
Temperature	Range & Sensor	Ambient + 3°C~80°C (0,1°C), RTD(Pt-100), PID control							
Pressure Control	The back pressure is controlled manually and/or automatically								
Control	Panel	17" Touch Screen Color LCD, Built-in industrial computer							
	Control program	Process, Trend, Calibration Mode control display, all the programmed setting values will be automatically stored into the built-in memory.							
	Processor	PLC system							
	Setting	Automated sterilization process Automated fermentation program Valve control							
	Communications	Memory stick, printer, computer, etc with the RS-232C communication protocol							
Data Storage	pH, DO, Temperature, Pressure and Agitation speed								

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- Vessel Ratio (H/D)
 - 2:1
 - Total Vol. L / Working Vol. L
- Vessel Pressure Certification Stamp
 - No Stamp / ASME Stamp
- Vessel Surface finish
 - Inside / outside : 0,4~0,8 μ m / 0,8~1,2 μ m
- Sight glass 1 ea : Oblong type
- Top Plate
 - Exhaust/Port / Pressure Gauge with diaphragm / Pressure Transmitter / Foam Sensor / Septum Port 4ea / Inoculum Port
- Upper Port
 - Aeration / Aeration by pass / Circle Light / Rupture Disk
- Lower Port
 - Sampling / pH / DO / Temp / Etc.
- Harvest

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




- According to customer requirements and the intended use by changing some standards it can be produced.
- The control software also can be produced according to customer requirements.

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Specification


Light intensity control range	10 ~ 100%
Light module	Fluorescent LED Lamp
Light Color	White, RED, Blue, Green
Max. Lux	Max. 10000lux
Sensor Control & Recording	pH, DO, Temp, Foam, Pressure etc.
Sparger & Gas Mode	Micro sparger, 2 gas and 4 gas mixer(optional)
Filter & Condenser	0.2µm air filter cartridge, SUS 316
Temp Control System	Cartridge heater or Chiller
Feeding Control System	Peristaltic pump Control Acid, Base, Foam, Feed
Light Control Mode	On/off or set of the lighting Controller
Power Voltage	110V~220V, 50/60Hz, 15A

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Cat. NO.	Description
FSBC-PP-PB01	Photo Bioreactor Autoclavable Fermentor 30L ~ 200L

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Fluorescent
LED lamp Case

BORULU SÜREKLİ SANTRİFÜJ

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Model	CT-Station 2	CT-Station 6	CT-Station 9
Special Feature	Open the front open (SUS Type)	Open the front open (SUS Type)	Open the front open (SUS Type)
Speed	21,000rpm	19,000rpm	15,000rpm
Centrifugal Effect(G)	22,500g	21,000g	17,000g
Throughput	500l /h	1,500l /h	3,000l /h
Bowl Capacity	2.7l	6l	9l
Weight	300Kg	480Kg	550Kg

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Ordering Information


Cat. NO.	Description
CC-CT-2	Industrial Continuous Centrifuge
CC-CT-6	Industrial Continuous Centrifuge
CC-CT-9	Industrial Continuous Centrifuge

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



- | | |
|----|----------------------|
| 1 | FEED IN |
| 2 | DRAG FIXED HANDLE |
| 3 | DRAG |
| 4 | FEED NOZZLE |
| 5 | BOWL |
| 6 | THREE WING |
| 7 | DISCHARGE PORT |
| 8 | COVER |
| 9 | SPINDLE COVER NUT |
| 10 | SPINDLE |
| 11 | BEARING |
| 12 | PULLEY |
| 13 | TENSION PULLEY |
| 14 | MOTOR |
| 15 | COOLING COIL OPTIONS |
| 16 | FRAME |

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DISK SÜREKLİ SANTRİFÜJ

CDI-İstasyon Seris

Model	Capacity(L/h)	Maximum Rotational Speed(rpm)	Dimension
CDI-STATION 15D	400~1000	7000	950x950x1250mm
CDI-STATION 35D	1000~3000	7000	1400x1400x1450mm
CDI-STATION 85D	3000~7000	7000	1750x1450x1850mm

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
Ordering Information

Cat. NO.	Description
CC-CDI-15D	Industrial Continuous Centrifuge
CC-CDI-35D	Industrial Continuous Centrifuge
CC-CDI-85D	Industrial Continuous Centrifuge

1. Feeding Water
2. Product inlet
3. Discharged clear liquid phase
4. Centripetal pump
5. Disc group
6. Discharged heavy liquid phase
7. Residue space
8. Sparation disc

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
CENTRION Co.,Ltd

ANA OFİS

Bucheon Techno Park Ssangyong 3cha 397,
Seokcheon-ro, Ojeong-gu, Bucheon-si, Gyeonggi-do, Rep. of KOREA

Tel +82 32-624-2060/2160

Fax +82 32-624-2059

Website www.centrion.co.kr

Email sales@centrion.co.kr / centrion@centrion.co.kr

Yaşam Bilimleri İle Temas Halinde