





biopexfill®

biopex Fill



biopexFill® A Possibility for Everyone's Beauty



biopexFill®

Established in 2020, Glopex Co., Ltd is fast growing in the global health and aesthetic markets, aiming to be a developer and manufacturer of biomedical and biopharmaceutical products. Glopex®'s main business area focuses on the biomedical field to improve overall health, life, and beauty.

With continuous in-house research and development, we strive to develop various advanced products, including HA fillers, skin boosters, botulinum toxins, and other aesthetic products with our exclusive technologies.

Glopex® has innovated its biopexFill® HA gel lines, featuring a variety of injectable cross-linked hyaluronic acid gels in 1mL and 10mL syringes, as well as 60mL vials with distinctive applications for the body. These products come in various formulations and densities, contributing to increased volume and fullness of the skin

biopexfill® is a distinctive performance hyaluronic acid gel, showing the most natural and satisfactory results.

Glopex® distinguishes itself with exceptional customer service, providing exclusive benefits and rewards to members of the biopexFill® club.

biopexFill® 1mL line

biopexFill® is Glopex®'s first line of HA gels consisting of biopexFill® FINE, biopexFill® DEEP, biopexFill® SUB-Q and biopexFill® ULTIMATE. We offer four distinct viscoelasticity products for customized wrinkle correction and contouring treatments.

biopexFill®		biopexFill® FINE FINE biopex Fill® , weare related to the second seco	biopexFill® DEEP biopexFill®,	biopexFill® SUB-Q SUB-Q biopexFill® >
	Dosage	1m.L	1mL	1mL
	HA Concebtration	24mg/mL	24mg/mL	24mg/mL
	Lidocaine Content	3mg/mL	3mg/mL	3mg/mL
	Cross-link Agent	BDDE	BDDE	BDDE
\	Needle Size	30G	27G	25G
	Depth of Injection	Superficial to mid dermis	Superficial to mid dermis	Deep dermis
	Indication	 Forehead fine lines Lips lines (Smoker's line) Horizontal neck lines Undereye Superficial fine lines Glabellar fine lines 	 Lips contour and Downturned lip corners Marionette lines Nasolabial fold Malar area Midface augmentation Deep frontal lines 	NoseNasolabial foldPrejowl sulcusDeep nasolabia linesTemporal
	Volume creating Capability	••000	•••00	••••

biopexFill® 10mL line

The biopexFill® second-line product is a 10mL HA gel with ideal viscoelasticity for convenient, versatile procedures. It's included connector allows easy division into lower-dosage syringes, minimizing contamination risks.


biopexFill® ULTIMATE	biopexFill® LIPS-S10	biopexFill® FACE-S10	biopexFill® BODY-V60
biopex Filli ,	195-510 Lipota da est atrada de la state	MCCSB biopex Fill bipodes for an including 185 mt	biopex fill
1mL	10mL	10mL	60mL
24mg/mL	24mg/mL	24mg/mL	20mg/mL
3mg/mL	3mg/mL	3mg/mL	Whitout Lidocaine
BDDE	BDDE	BDDE	BDDE
23G			
Supraperiosteal	Mid to deep dermis	Supraperiosteal	Deep dermis to subcutaneus
- Chin augmentation	- Cheek augmentation	- Facial contours	- Hand rejuvention
- Jaw line	- Lips line, contour,	- Cheek augmentation	- Breast and buttock
- Angulus mandibula	and downturned lip	- Nose	augmentation
	corners	- Nasolabial fold	- Hip dips correction
	- Temples	- Marionette lines	- Cellulite reduction
			- Calves augementation
			- Male genital
•••••	••000	•••••	•••00

biopexFill® 60mL


The "BODY-V60" is a unique 60mL HA gel designed to enhance the desirable appearance of breasts, legs, and buttocks. Lidocaine-free, it is suitable for body injections requiring larger volumes of HA gel.

High-Density Hyaluronic Acid Gel

biopexFill® is manufactured of high-density Hyaluronic Acid gel using unique cross-linking technology to create a higher viscoelasticity gel that is safe to use.




Advantage of biopexFill® products rheological versatility




biopexFill® products provide a various selection of hyaluronic acid fillers, ensuring predictable and natural results to meet the individual

biopexFill® Safety & Purification



Modification Degree (MoD)

A lower Modification Degree makes it more similar to the natural properties of hyaluronic acid, ensuring a safer gel. biopexFill® has a low MoD, providing safety with confidence.



• Removal of Residual BDDE

biopexFill® refinement method minimizes the residual BDDE through a two-week purification process to improve the product's safety.



1.0 mL Line FINE / DEEP / SUB-Q / ULTIMATE

biopexFill® FINE offers the smoothest viscoelasticity, ideal for fine wrinkle treatments like forehead lines and neck wrinkles, as well as soft lip-line treatments.

biopexFill® DEEP is suitable for immediate correction of deep wrinkles, marionette lines, cheeks, chin, temple, lip contour, and scar treatment.

biopexFill® SUB-Q, with high gel density, is best for restoring fullness in nasolabial folds, jawline, facial contour, and nose, enhancing facial profile shape and definition.

biopexFill® ULTIMATE, the firmest with the highest viscoelasticity, provides immediate effects for facial contour correction, volumization, and is designed for treatments like chin augmentation and jawline enhancement.

biopex Fill®





10mL Line LIPS-S10 / FACE-S10

biopexFill® LIPS-S10 is a product specialized for lip treatments with a soft viscoelasticity. It is suitable not only for lip line treatments and lip augmentation but also for cheek augmentation and temple.

biopexFill® FACE-S10 is optimized for facial contouring and augmentation, featuring a high-viscoelasticity gel. This characteristic makes it effective for filling wrinkles such as nasolabial fold, marionette lines, and nose shaping.





60mL Line BODY-V60

biopexFill® BODY-V60 can be the perfect option for restoring volume in areas of significant volume loss, such as the breasts, lower legs, buttocks, hip dips, and genital areas. Additionally, it is an effective product for hand rejuvenation and cellulite reduction.







A Possibility for Everyone's Beauty



Glopex Co., Ltd.

Headquarters: A2801, R&D: B2065

 $IX\ Tower, Dongtan cheom dans an eop\ 1-ro\ 27, Hwaseong-si, Gyeonggi-do, 18469, Republic\ of\ Korea$

Factory: 27-10, Gyeonggi-daero 658beon-gil, Osan-si, Gyeonggi-do, 18112, Republic of Korea

Websites: www.glopex.kr / www.biopexfill.com